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15. Referenzen

1. American Cancer Society. www.cancer.org, accessed March 2006.
2. Fleshner NE, HW Herr, AK Stewart, GP Murphy, C Mettlin, HR Menck. The National Cancer Data Base Report on Bladder Carcinoma Cancer. 78: 1505-1513, 1996.

Rev. 3, 2025-11-07 OM

1. Intended Use

The NADAL® Nuclear Matrix Protein 22 Test is an *in-vitro* immunoassay intended for the qualitative detection of the nuclear mitotic apparatus protein (NuMA) - an abundant component of nuclear matrix proteins - in the urine of individuals with risk factors or symptoms of bladder cancer, or with a history of bladder cancer. This assay provides a preliminary result only. Clinical expertise and professional judgment should be sought to further evaluate the result of the test.

2. Introduction and Clinical Significance

It was estimated by the American Cancer Society that there would be 61,420 new cases of bladder cancer in the United States in 2006, with approximately 13,060 of these patients dying from the disease.¹ The majority of bladder cancer cases occur in individuals over 50 years of age, and over 70% of these cases occur in men. The early diagnosis of cancer has a dramatic impact on prognosis: the five-year survival rate for patients with bladder cancer is 94% - if the tumor is diagnosed at a superficial stage; cancers detected when they have become invasive or metastatic have five-year survival rates of 49% and 6% respectively. 75% to 80% of bladder tumors present as superficial disease (Ta, Tis, T1), and the majority of patients with superficial tumors (50-80%) will experience recurrence within six to twelve months. Approximately 10-20% of those tumors will progress in stage and grade. The predominant risk factors for recurrence are large tumor size and high tumor grade. Most commonly, recurrence is at the site of the original tumor resection, but new tumor formations can also develop at other sites within the bladder.²

3. Test Principle

The NADAL® Nuclear Matrix Protein 22 Test is for the detection of NMP22 through the visual interpretation of colour development on the internal test strip. Anti-NMP22 antibodies are immobilised in the test line region (T) of the membrane. During the test, the specimen reacts with anti-NMP22 antibodies which are conjugated to coloured particles and precoated onto the conjugate pad of the test cassette. The mixture then migrates along the membrane by capillary action and interacts with the reagents on the membrane. If there is a sufficient amount of NMP22 in the specimen, a coloured line will develop in the test line region (T) of the membrane. The presence of this coloured line indicates a positive result, while its absence indicates a negative result. The formation of a coloured line in the control line region (C) serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred. If the control line does not appear, the test result is to be considered invalid.

4. Reagents and Materials Supplied

- 2/10/20 NADAL® Nuclear Matrix Protein 22 test cassettes
- 2/10/20 disposable pipettes
- 1 package insert

The NADAL® Nuclear Matrix Protein 22 Test contains anti-NMP22 antibodies conjugated to coloured particles and anti-NMP22 antibodies precoated onto the membrane.

5. Additional Materials Required

- Timer

6. Storage and Stability

Test kits should be stored at room temperature or refrigerated (2-30°C). Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze test kits. Do not use tests beyond the expiry date indicated on the packaging.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Do not use the test beyond the expiry date indicated on the packaging. Do not use the test if the foil pouch is damaged. Do not reuse tests.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Carefully read through the entire instructions for use prior to testing.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- Humidity or temperature can adversely affect test results.
- Used testing materials should be disposed of in accordance with local regulations.

8. Specimen Collection and Preparation

Voided urine is required for use with the NADAL® Nuclear Matrix Protein 22 Test. Barbotage specimens, catheterised urine specimens, serum, plasma or whole blood should not be used. Urine should be collected without preservatives or fixatives in a clean urine collection cup and labeled appropriately. If the urine specimen is to be used for other, additional tests, remove an aliquot of the specimen (a minimum of 2 mL) for this test to avoid contamination. Urine specimens can be stored at room temperature for up to 2 hours after collection. The test should not be used on individuals with indwelling urinary tract devices (such as stents) or those who have had a total cystectomy. Use only plastic containers for urine specimen collection. Do not use paper or foam cups. Urine collection during active chemo-, immuno-, or radiation therapy may cause positive test results. Some patients with active prostate cancer may yield positive results with the NADAL® Nuclear Matrix Protein 22 Test. Before carrying out the tests, the physician should allow ample time for recovery from trauma to the bladder or urinary tract resulting from surgery, biopsies, etc.

9. Test Procedure

Bring tests, specimens and/or controls to room temperature (15-30°C) prior to testing.

1. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed within one hour.
2. Place the test cassette on a clean and level surface.
3. Holding the pipette vertically, transfer 2-3 drops of the urine specimen to the specimen well (S) of the test cassette.

Avoid trapping air bubbles in the specimen well (S).

4. Start the timer.
5. Wait for the coloured line(s) to appear. Read the result in 10-15 minutes. Do not interpret results before 10 minutes or after 15 minutes. **It is important to check that the background is clear before reading the result.**

- Results of the NADAL® Nuclear Matrix Protein 22 Test should not be interpreted as evidence of the presence or absence of a malignant disease in the bladder without confirmation using other diagnostic procedures. To ensure good patient management, other clinically accepted tests and procedures should be considered before diagnosing a disease.

13. Performance Characteristics

Analytical Sensitivity

The NADAL® Nuclear Matrix Protein 22 Test has a sensitivity of 20 U/mL urine.

Diagnostic sensitivity and specificity

To obtain diagnostic sensitivity and specificity data, urine from 32 patients with biopsy-confirmed bladder cancer and 112 healthy individuals with no bladder tumors was tested using the NADAL® Nuclear Matrix Protein 22 Test.

The results are presented in the following table:

NADAL® Nuclear Matrix Protein 22 Test		Biopsy		
		Positive	Negative	Total
	Positive	30	9	39
	Negative	2	103	105
	Total	32	112	144

Diagnostic sensitivity: 93.8% (30/32) (95% CI, 77.78% - 98.91%)

Diagnostic specificity: 92.0% (103/112) (95% CI, 84.88% - 96.03%)

Interference study

Potentially interfering substances were spiked with negative and positive urine specimens with concentrations listed in the table below. None of the tested substances showed any interference with the NADAL® Nuclear Matrix Protein 22 Test.

Substance	Concentration	Substance	Concentration
Control	N/A	Trimethoprim	0.5 mg/mL
Whole blood	0.6%(v/v)	Flutamide	0.8 mg/mL
Protein	0.2 mg/mL	Bilirubin	0.008 mg/mL
Glucose	0.1 mg/mL	Nitrofurantoin	0.5 mg/mL
Mitomycin C	0.2 mg/mL	IgG	0.1 mg/mL
Acetaminophen	0.2 mg/mL	Uric acid	2 mg/mL
Sodium salicylate	0.2 mg/mL	Haemoglobin	1 mg/mL
Aspirin	0.2 mg/mL	Ciprofloxacin	1 mg/mL
Ibuprofen	0.2 mg/mL	Levofloxacin	1.5 mg/mL
Ampicillin	0.2 mg/mL	Isovue®	0.175 mg/mL
Vitamin C	0.2 mg/mL	Thiotepa	0.6 mg/mL
Sodium chloride	3 mg/mL	Finasteride	0.025 mg/mL
Ethanol	1%(v/v)	Phenazopyridine hydrochloride	0.7 mg/mL
Tetracycline	0.2 mg/mL	Doxorubicin hydrochloride	0.1 mg/mL

14. Serious incident reporting

In the event of any serious incidents related to the performance of the NADAL® Nuclear Matrix Protein 22 Test, please inform nal von minden GmbH and the competent

10. Result Interpretation

Positive:

A coloured line develops in the control line region (C) and another coloured line develops in the test line region (T).

Negative:

A coloured line develops in the control line region (C). No line develops in the test line region (T).

Invalid:

The control line (C) fails to appear. Insufficient specimen volume or incorrect operating procedure are the most likely reasons for control line failure. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

11. Quality Control

An internal procedural control is included in the test cassette: A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this test kit. However, it is recommended that positive and negative controls be obtained from a local competent authority and tested as *good laboratory practice* to confirm the test procedure and to verify proper test performance.

12. Limitations

- Results of the NADAL® Nuclear Matrix Protein 22 Test should not be used as absolute evidence of the presence or absence of bladder cancer. Any disease which may trigger the presence of nuclear matrix proteins in urine may cause a positive test result. Positive results have been observed in some patients with benign urological diseases, prostate cancer and those undergoing active cancer treatment.
- Results of the NADAL® Nuclear Matrix Protein 22 Test should be interpreted only in conjunction with information obtained from the clinical evaluation of the patient and other diagnostic procedures.

NADAL® Nuclear Matrix Protein 22 Test

(Ref. 562011N-02/ 562011N-10/ 562011N-20)

authority immediately. If still possible, do **not** dispose of the used test and the corresponding test kit components.

15. References

1. American Cancer Society. www.cancer.org, accessed March 2006.
2. Fleshner NE, HW Herr, AK Stewart, GP Murphy, C Mettlin, HR Menck. The National Cancer Data Base Report on Bladder Carcinoma Cancer. 78: 1505-1513, 1996.

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Zhejiang Orient Gene Biotech Co., Ltd.