natural airway management


Instructions For Use

www.i-gel.com

STERILE EO

Rx ONLY

natural airway management

i-gel® supraglottic airway	
en	Instructions for use
fr	Mode d'emploi
de	Gebrauchsanweisung
es	Instrucciones de uso
pt	Instruções de utilização
it	Istruzioni per l'uso
nl	Gebruiksaanwijzing
no	Instruksjoner for bruk
fi	Käyttäjän opas
sv	Bruksanvisning
da	Brugsanvisning
el	Οδηγίες χρήσης
lt	Naudojimo vadovas
pl	Instrukcja użytkowania
ru	Инструкция по применению
cs	Návod k použití
hu	Használati utasítás
sl	Navodila za uporabo
lv	Lietošanas instrukcija
et	Kasutamisujuhised
bg	Инструкции за употреба
ro	Instrucțiuni de utilizare
hr	Upute za upotrebu
tr	Kullanma kılavuzu
ja	取扱説明書
zh-s	使用说明
zh-t	使用說明
ar	تعليمات الاستخدام

Figure 1


Figure 2



Figure 3

The medical device is not supplied individually. Only one copy of the instructions for use is provided per box and should therefore be kept in an accessible location for all users.

NOTE: Distribute instructions to all product locations and users.

These instructions contain important information for the safe use of the product. Read these instructions for use in their entirety, including warnings and cautions before using this product.

WARNING: A WARNING statement provides important information about a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation, which if not avoided could result in minor or moderate injury.

Product name & description

i-gel® supraglottic airway (paediatric sizes). Product codes: 8225000, 8202000, 8215000 and 8201000

Key components (Figure 2)

- 1.1 Distal end of gastric channel (8225000, 8202000 and 8215000)
- 1.2 Epiglottic rest
- 1.3 Proximal end of gastric channel (8225000, 8202000 and 8215000)
- 1.4 Soft non-inflatable cuff
- 1.5 Buccal cavity stabiliser
- 1.6 Integral bite block
- 1.7 15mm connector

Product positioning (Figure 3)

Key:

- 2.1 Oesophagus
- 2.2 Trachea
- 2.3 Vocal cords/folds
- 2.4 Thyroid cartilage
- 2.5 Laryngeal Inlet
- 2.6 Sealing mechanism. Non-inflatable cuff
- 2.7 Ventilatory opening
- 2.8 Hyoid bone
- 2.9 Epiglottis
- 2.10 Tongue
- 2.11 Nominal depth of insertion mark
- 2.12 Incisor
- 2.13 Nominal length of the internal pathway (refer to technical data chart below)

Intended use

Securing and maintaining a patent airway. To provide a conduit for the delivery of respiratory gases to the patient.

Indications

The paediatric i-gel® is indicated in securing and maintaining a patent airway in routine and emergency anaesthetics for operations of fasted patients during spontaneous or intermittent positive pressure ventilation (IPPV) and during resuscitation of the unconscious patient, by personnel who are suitably trained and experienced in the use of airway management techniques and devices.

Contraindications

- Non-fasted patients for routine and emergency anaesthetic procedures.
- Trismus, limited mouth opening, pharyngo-perilaryngeal abscess, trauma or mass.
- Do not allow peak airway pressure of ventilation to exceed 40cm H₂O.
- Do not use excessive force to insert the device or gastric tube.
- Inadequate levels of anaesthesia which may lead to coughing, bucking, excessive salivation, retching, laryngospasm or breath holding thus complicating the anaesthetic outcome.
- Do not leave the device in situ for more than 4 hours.
- Do not reuse or attempt to reprocess the i-gel®.
- Patients with any condition which may increase the risk of a full stomach e.g. hiatus hernia, sepsis, morbid obesity, pregnancy or a history of upper gastro-intestinal surgery etc.
- Use on a conscious/semi-conscious patient in an emergency setting.

Patient Category

The product range is suitable for use with all patient populations. The product range is suitable for use with paediatric patients with a body weight as indicated on the individual product pack and the connector of the device. Please refer to the technical data table for size selection. Whilst size selection on a weight basis should be applicable to the majority of patients, individual anatomical variations mean the weight guidance provided should always be considered in conjunction with a clinical assessment of the patient's anatomy.

Warnings

- The i-gel® is supplied in a cage pack to ensure the device is retained in the correct flexion prior to use and also acts as a base for lubrication. The i-gel® must always be separated from the cage pack prior to insertion. The cage is not an introducer and must never be inserted into the patient's mouth.
- i-gel® must be lubricated according to the instructions for use.
- Remove dentures or removable plates from the mouth before attempting insertion of the device.
- The patient should always be in the 'sniffing-the-morning air' position prior to insertion with the assistant helping to open the patient's mouth, unless head/neck movements are considered inadvisable or are contraindicated.
- Optimum depth of anaesthetic must be achieved prior to attempting insertion (i.e. absence of eyelash reflex, easy up and down movement of the lower jaw, no reaction to pressure applied to both angles of the mandible.)
- The leading edge of the i-gel® tip must follow the curvature of the patient's hard palate upon insertion.
- Do not apply excessive force on the device during insertion. It is not necessary to insert fingers or thumbs into the patient's mouth during the process of inserting the device.
- If there is a failure to achieve complete insertion after utilising the standard insertion technique and a jaw thrust, deep rotation or triple manoeuvre has also failed, then the device should be inserted under direct vision by laryngoscopy or one size smaller device should be used.
- After insertion, i-gel® should be taped from maxilla-to-maxilla. The i-gel® should be tied as well as taped during surgical procedures that require the patient to be in a lateral or prone position. The i-gel® may also be tied in place in circumstances where tape is considered inadequate or ineffective.
- Excessive air leak during manual ventilation is primarily due to either sub-optimal depth of anaesthesia or sub-optimal depth of i-gel® insertion.
- Particular care should be taken with patients who have an ASA or Mallampati score of III and above or who have fragile and vulnerable dental work, in accordance with recognised airway management practices and techniques.
- As with all supraglottic airways, it is important to ensure the correct size of device is used, lubrication is optimal, the device is inserted and positioned correctly and regularly checked intraoperatively in order to reduce the potential for nerve damage, tongue numbness, cyanosis and other potential complications.
- As with any airway device, when using the i-gel®, it is essential to ensure there is always an adequate expiratory pathway for gas to avoid the potential for barotrauma.
- A gastric tube should not be passed down the gastric channel of the i-gel® if there is an excessive leak through the gastric channel, or in cases of oesophageal trauma, or if there are oesophageal varices or evidence of upper gastro-intestinal bleed, there is a history of upper gastro-intestinal surgery or the patient has bleeding/clotting abnormalities.
- Insertion of a gastric tube in the presence of inadequate levels of anaesthesia can lead to coughing, bucking, excessive salivation, laryngospasm or breath holding.
- Do not attempt to forcibly remove the device if the patient is biting on it. Wait until the patient, on vocal command, has fully opened their mouth or opens their mouth spontaneously.
- Supraglottic airways may be flammable in the presence of lasers and electrocautery equipment.
- The i-gel® incorporates a gastric channel (except 8201000 - Size 1), which, when the device is correctly placed, may provide some protection against aspiration of gastric contents. However, the level of protection will vary and cannot be considered definitive.
- In order to reduce the risk of fire, it is essential the user fully complies with the guidance provided by the European Resuscitation Council (ERC), the American Heart Association (AHA) or other appropriate local resuscitation council regarding the safe use of oxygen during defibrillation.
- Whilst size selection on a weight basis should be applicable to the majority of patients, individual anatomical variations mean the weight guidance provided should always be considered in conjunction with a clinical assessment of the patient's anatomy. Patients with cylindrical necks or wide thyroid/cricoid cartilages may require a larger size i-gel® than would normally be recommended on a weight basis. Equally, patients with a broad or stocky neck or smaller thyroid/cricoid cartilage, may require a smaller size i-gel® than would normally be recommended on a weight basis. Patients with central obesity, where the main weight distribution is around the abdomen and hips, might in practice require an i-gel® of a size commensurate with the ideal body weight for their height rather than their actual body weight.
- In patients with extremely low lung compliance, where it is not possible to deliver adequate ventilation with i-gel®, an alternative airway device should be used.
- For patients undergoing Electroconvulsive Therapy (ECT) treatment, a separate bite block specifically designed to be used during ECT must always be used in conjunction with the i-gel®.

Cautions

- Always perform the pre-use checks included in these instructions for use.
- The patency of the device should be checked regularly and particularly after any change in the patients head or neck position.
- i-gel® is a single use device. Do not reuse. Reuse poses a threat to patient safety in terms of cross infection, device malfunction and breakage. Cleaning can lead to device malfunction and breakage.
- The size 1 i-gel® does not have a gastric channel.
- The i-gel® should always be used in accordance with recognised airway management practice for supraglottic airway devices.

Size Selection

Select the appropriate size i-gel® by assessing the patient's anatomy and referring to the patient size and patient weight guidance in the technical data chart below. The i-gel® cuff may look smaller than traditional supraglottic devices with an inflatable cuff of the same numerical size.

Pre-use checks

- Inspect the packaging and ensure it is not damaged prior to opening.
- Inspect the device carefully, check the airway is patent and confirm there are no foreign bodies or a BOLUS of lubricant obstructing the distal opening of the airway or gastric channel.
- Carefully inspect inside the bowl of the device, ensuring surfaces are smooth and intact and also that the gastric channel is patent.
- Discard the device if the airway tube or the body of the device looks abnormal or deformed.
- Check the 15mm connector fits the patient connection.

Pre-insertion preparation

1. Always wear gloves
2. Open the i-gel® package, and take out the cage pack containing the device (Figure 1, picture 1)
3. In the final minute of pre-oxygenation, open the cage pack and transfer the device into the lid of the cage. Place a small bolus of a water based lubricant onto the smooth inner surface ready for use. Do not use silicone based lubricants. (Figure 1, picture 2,3 and 4)
4. Grasp the i-gel® along the integral bite block and lubricate the back, sides and front of the cuff with a thin layer of lubricant. Ensure after lubrication has been completed, that no BOLUS of lubricant remains in the bowl of the cuff or elsewhere on the device. Avoid touching the cuff of the device with your hands (Figure 1, picture 5, 6, 7 and 8).
5. Place the i-gel® back into the cage pack in preparation for insertion (Figure 1, picture 9).

Recommended insertion technique

A proficient user can achieve insertion of the i-gel® in less than 5 seconds.

1. Remove the i-gel® from the cage pack and grasp the lubricated i-gel® firmly along the integral bite block. Position the device so that the i-gel® cuff outlet is facing towards the chin of the patient. (Figure 1, picture 10)
2. The patient should be in the 'sniffing the morning air' position (Figure 1, picture 10) with head extended and neck flexed. The chin should be gently pressed down by an assistant before proceeding to insert the i-gel®.
3. Introduce the leading soft tip into the mouth of the patient in a direction towards the hard palate.
4. Glide the device downwards and backwards along the hard palate with a continuous but gentle push until a definitive resistance is felt.
5. At this point the tip of the airway should be located into the upper oesophageal opening (Figure 1, picture 11a) and the cuff should be located against the laryngeal framework (Figure 1, picture 11b). The incisors should be resting on the integral bite-block (Figure 1, picture 11c).
6. i-gel® should be taped down from 'maxilla to maxilla' (Figure 1, picture 12).
7. If required, an appropriate size gastric tube may be passed down the gastric channel (except size 1.0). Please see the technical data chart below for guidance.

Removal

- Once consciousness is regained and protective reflexes such as coughing and swallowing have returned, gently suction around the airway device in the pharynx and hypopharynx. Once the patient is awake or easily arousable with vocal commands, the i-gel® can safely be removed by asking the patient to open his/her mouth.
- In patients with the possibility of a heightened gag reflex (i.e. smokers, asthmatics or patients with COPD), i-gel® should be removed in deeper planes of anaesthesia and, after removal, the airway maintained with a Guedel airway until protective reflexes have returned and the patient becomes arousable.

Storage conditions

Store in a dry place out of direct sunlight.

Disposal

Following use, the product must be disposed of in accordance with local hospital infection control and waste disposal regulations.

Technical data

Product code	Size	Patient weight (kg)	Internal volume (ml)	Nominal length of internal pathway (mm)	Maximum size of ET tube that can be inserted down the airway channel (mm)	Maximum size of gastric tube that can be passed down the gastric channel (FG)	Minimum interdental gap for insertion (mm)	Pressure drop
8201000	1	2 - 5	5.5	119	3.0	Not applicable	13.9	< 0.6cmH ₂ O @ 15L/min
8215000	1.5	5 - 12	8.0	145	4.0	10	17.6	< 0.3cmH ₂ O @ 30L/min
8202000	2	10 - 25	13.9	160	5.0	12	21.3	< 0.3cmH ₂ O @ 30L/min
8225000	2.5	25 - 35	20.7	182	5.0	12	24.9	< 0.5cmH ₂ O @ 60L/min

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Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK