ENGLISH

INSTRUCTIONS FOR USE

Follow instructions for use, manufacturer guidelines, and institution procedures for flush administration.

DESCRIPTION

BD PosiFlush™ XS Syringe is a ready to use sterile medical device (according to regulation (eu) 2017/745 of the European Parliament and of the Council). It is a polypropylene syringe containing sterile and non-pyrogenic isotonic 0.9% sodium chloride solution. The contents of our unopened or undamaged blister packages are guaranteed to be sterile, non-toxic, and non-pyrogenic.

INTENDED USE/INDICATIONS FOR USE

- BD PosiFlush™ XS Syringes are intended to be used FOR FLUSHING ONLY of in-situ peripheral intravenous catheters (PIVCs), peripherally inserted central catheters (PICCs), central venous catheters (CVCs), and implanted venous access ports.
- BD PosiFlush™ XS Syringe is not intended for dry product reconstitution, for medication administration, or when intravenous therapy with sodium chloride is indicated.
- Using aseptic technique, **BD PosiFlush™ XS Syringe can be used on a sterile field.**

INTENDED PATIENT POPULATION

BD PosiFlush™ XS Syringe to be used with patients with in-situ peripheral intravenous catheters (PIVCs), peripherally inserted central catheters (PICCs), central venous catheters (CVCs), and implanted venous access ports.

INTENDED USER

BD PosiFlush™ XS Syringe is to be used by healthcare professionals experienced in vascular access and the use of these devices.

The BD PosiFlush™ XS Syringes are manual devices that may be operated by a large range of people with various human characteristics, including gender, age, and weight.

CONTRAINDICATIONS

- Do not use in patients suffering from hypernatremia and fluid retention when the administration of sodium or chloride could be detrimental.
- Do not use BD PosiFlush™ XS Syringe if a patient has is known allergy to any of its components, mainly to the 0.9% sodium chloride solution, which may lead to an allergenic response resulting in anaphylaxis.

PERFORMANCE CHARACTERISTICS

- Interoperability with Vascular Access Devices (VAD).
- All syringes of 3 mL, 5 mL, and 10 mL have a maximum 10 mL diameter syringe barrel and therefore the flush pressure is equivalent for all sizes.

WARNINGS

- Do not use if unit package or content is damaged.
- Do not use if product is past the labeled expiration temperature.
- Verify the expiration date on the product package or label. Do not use if product has expired.
- Do not use if syringe tip cap or stopper is damaged.
- Do not use if solution is cloudy, contains a precipitate, or has any type of suspended particulate matter.
- Do not remove the luer stopper due to risk of infection or other illness/injury.
- Sodium chloride has a potential choking hazard. After use, discard small parts according to your facility protocol.
- Do not re-sterilize.
- Possible complications and/or adverse reactions associated with flushing may include sepsis, infection (localized/systemic), mucocutaneous blood exposure, exposure to bloodborne pathogens, air embolism, particulate embolism from clots, catheter occlusion, and damage to the vessel wall. Single blister pack consists of paper and Polypropylene – Polyamide film – Polypropylene (PPE).

PRECAUTIONS

- Store at controlled temperature (15-25°C). Excursions permitted to 30°C. Do not leave in temperature extremes.
- Check with drug manufacturer instructions for use to ensure compatibility before flushing. Do not use BD PosiFlush™ XS if 0.9% sodium chloride solution is not compatible; follow the drug manufacturer instructions for flushing practices, or first flush the vascular access device (VAD) with a compatible solution such as 5% dextrose or water to remove traces of the medication in accordance with manufacturer and institution policies.
- Use aseptic technique in accordance to specific medical conditions, treatment needs, age, and weight (that may require restricted sodium or fluid intake when deciding to flush with 0.9% sodium chloride injection. When flushes should be limited due to the patient's condition, do not exceed fluid intake guidelines.
- Use an appropriate volume and technique of flushing with BD PosiFlush™ XS Syringe in pregnant and lactating women.
- BD PosiFlush™ XS Syringe is not to be used with ISO luer compliant components for intravenous applications.
- Do not use if syringe is damaged or already used product.
- EU Only: Users should report any serious incident related to the device to the BD manufacturer and competent authority.

DIRECTIONS FOR USE

Single use, single patient device. To ensure safe medication preparation and administration, clinicians should practice the "7 rights" of medication administration: right patient, right medication, right dose, right route, right reason and right documentation.

Use aseptic technique and precautions to flush before and after each medication, fluid administration, or blood sampling and at regular intervals during continuous infusions to maintain the patency of vascular access device, manufacturer, and institution policies.

Follow disinfection protocols as per institution and manufacturer guidelines.

Use aseptic technique throughout the procedure.

1. Open blister pack and remove syringe using sterile technique, place syringe on a sterile field if applicable.
2. Check that syringe is not damaged. Inspect clarity of solution. (Fig. 1a)
3. Depress plunger with tip cap on to release the stopper seal. (Fig. 1b)
4. Unscrew tip cap from the syringe ensuring that there is no touch contamination of the connection. (Fig. 2)
5. Push syringe plunger to expel the air. (Fig. 3)
6. Connect BD PosiFlush™ XS syringe to vascular access device, taking care that there is no touch contamination of the connection. (Fig. 4) Ensure product has expired.
7. Push syringe plunger to inject the required volume of saline following manufacturer's IFU. (Fig. 5) Inject at a rate that will result in a flush overpressure. If vascular device or drug resistance, it is recommended that a smaller syringe be used.
8. After use, dispose of in accordance with recognized procedure in your institution. (Fig. 6)

COMPOSITION PER UNIT

Polypropylene syringe with BD Luer-Lok™ tip.

BD PosiFlush™ XS Syringe is not formulated with Natural Rubber Latex. Saline solution: sodium chloride 9g/l (NaCl 0.9%), distilled sterile water to volume, preservative free.

Single blister pack consists of paper and Polypropylene – Polyamide film – Polypropylene (PPE).

BENEFICIAL CLINICAL

- The BD PosiFlush™ XS syringe is a pre-filled, single use 0.9% sodium chloride syringe that helps to improve clinical flushing and administering steps and time involved in the manual preparation of saline syringes.
- The BD PosiFlush™ XS syringe can help reduce risk of cross-contamination that may occur during manual preparation of saline flush syringes.

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