

Model Comparison



**SnapIT™
Personal**



Anodised, Single
Piece of Aluminium



Latex Free



Autoclavable



Re-Usable



Keyring for
On-the-go Use



**SnapIT™
Lite**



Lightweight
Plastic



Latex Free



Anti-Roll
Technology



Re-Usable



Keyring for
On-the-go Use



**SnapIT™
Trolley**



Anodised, Single
Piece of Aluminium



Latex Free



Autoclavable



Anti-Roll
Technology



Re-Usable



Keyring for
On-the-go Use



Sizing Guide

1-2ml

5-10ml

10-15ml

20-25ml

25-30ml

Regular:



Large:



Extra Large:



Ampoule sizing by millilitre is approximate



Sizing Guide

Regular Size


D = 7.15mm - 9.75mm
(0.2814" - 0.3838")

Large Size

D = 9.70mm - 12.45mm
(0.3818" - 0.4901")

Extra-Large Size

D = 12.40mm - 15.00mm
(0.4881" - 0.5905")



Maximum

L ≤ 35.00mm
(1.3779")